


EndoArt[®]

Artificial Endothelial Layer

Instructions for Use

 **EyeYon**
MEDICAL

INSTRUCTIONS FOR USE

EndoArt® Corneal Artificial Endothelial Layer

Cat. No.: EARRXXYYZZZ

1. IMPORTANT

Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.

The device is provided STERILE. DO NOT RESTERILIZE.

DO NOT USE the device if the package is open or damaged, or if the expiration date has passed.

The device is intended for SINGLE-USE only.

2. DEVICE DESCRIPTION

EndoArt® (Artificial Endothelial layer) is a permanent synthetic implant. The device serves as an inert physical barrier, blocking the ingress of fluid from the Anterior Chamber (AC) into the cornea. By sealing part of the inner cornea, the EndoArt® relieves corneal edema without functioning endothelium. The EndoArt® device is a copolymer of hydroxyethyl methacrylate and methyl methacrylate. It is dome shaped, optically clear, transparent, foldable implant, made of hydrophilic material designed to adhere to the posterior portion of the cornea.

EndoArt® specifications:

Diameter (DIA)	6.0 - 8.0 mm
Curvature Radius (R)	6.0 - 8.0 mm
Thickness (T)	0.03 - 0.05 mm
Power (D)	0.00 ± 0.25 D

3. INTENDED USE

The EndoArt® device is intended to be used as an endothelial prosthesis in patients with chronic corneal edema due to endothelial dysfunction.

4. INTENDED POPULATION

Patient aged ≥ 40 years with chronic corneal edema due to endothelial dysfunction.

5. CONTRAINDICATIONS

DO NOT USE the EndoArt® when any of the following conditions exist:

- Severely scarred cornea unfit for regular endothelial keratoplasty
- Band keratopathy
- History of corneal refractive surgery: LASIK, PRK, SMILE, Epi-keratoplasty (or any anterior lamellar keratoplasty) and Radial Keratotomy (RK)
- History of neurotrophic cornea
- History of persistent/ recurrent corneal erosion or persistent epithelial defect with difficulties in epithelial re-growth (erosion more than 2 months)
- History of posterior vitrectomy
- Active infection of the eye
- Diagnosis of severe dry eye such as Sjogren's Syndrome, GVHD
- Severe corneal distortion e.g.: keratoconus, pellucid marginal degeneration or keratoglobus (not relevant if the patient already underwent PKP)
- PUK (peripheral ulcerative keratitis)
- Phthisis bulbi or subject is at risk of developing phthisis.
- Irido Corneal Endothelial (ICE) Syndrome or any rare disease/ syndrome creating anatomical or physiological anomalies of the anterior chamber (e.g., corectopia)
- Subject with medically uncontrolled high intra ocular pressure
- Aphakia
- Severe Pseudophakodonesis / dislocation or partial dislocation of the IOL
- A large iris defect which could compromise intraoperative air bubble formation- aniridia more than 30%
- Recurrent posterior, intermediate or anterior uveitis.
- Limbal stem cell deficiency
- Severe ptosis
- Megalo cornea or micro cornea ($WTW \geq 13\text{mm}$ or $WTW \leq 10\text{mm}$)
- History of spontaneous corneal perforations

6. WARNINGS

- Do not use the device if the expiration date specified on the label has passed.
- Do not use the device if its sterility appears to be compromised.
- Single use only. Do not re-use, re-process or re-sterilize, all of which may compromise the structural integrity of the device and/or lead to device failure, which may result in patient injury, illness or death. Re-use, re-processing or re-sterilization also bears risk of contamination to the device and/or patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another.
- Eye rubbing after implantation may lead to implant detachment. Guide the patient to avoid touching/rubbing the eye during the first three months after implantation. Transparent eye shield wear is recommended during this period.
- Device removal may be required due to insufficient performance or adverse reaction.

7. PRECAUTIONS

- Carefully inspect the device prior to use, to verify that it has not been damaged and that its size, shape and condition are suitable for the intended procedure. DO NOT use if product appears to be damaged.
- DO NOT use the device if the package has been opened, as sterility may be compromised.
- Only physicians trained in performing corneal surgical ophthalmic procedures can use EndoArt®, after undergoing learning the methodology of implanting EndoArt® and obtaining certification.
- EndoArt® implantation should be in an environment appropriate for ocular surgery.

8. POTENTIAL ADVERSE REACTIONS

Complications associated with EndoArt® implantation include, but are not limited to:

Cornea

- Excessive corneal thinning
- Corneal melting or perforation
- Cornea abrasion, opacity, or haze
- Device detachment
- Worsening of corneal edema
- Epithelial thickening
- Hyperopic shift
- Vision deterioration
- Corneal neovascularization
- Bullae formation

Lens

- Cataract induction
- Implant's adhesion to hydrophilic IOL
- IOL dislocation

Glaucoma

- Increased ocular pressure/glaucoma
- Pupillary block

Dry eye

- Epiphora
- Eye burning
- Allergy

Other

- Infection/endophthalmitis
- Inflammation
- Retinal detachment
- Anterior or posterior synechiae
- Hypotony
- Hypoxia due to post op mandatory supine position in patients at risk (e.g. heart failure, COPD, obese)
- Orthostatic hypotension

9. INSTRUCTIONS FOR USE

Implantation of the EndoArt® is according to clinical discretion and common practice of the physician, and in compliance with the following instructions:

1. Prepare the patient according to standard surgical practice.
2. Position the patient under the surgical microscope.
3. Remove the epithelial layer for better visualization and corneal dehydration.
4. Create a primary entry (e.g. 1.2–3.0 mm) with keratome side port and secondary entry point (one or two as necessary) using a stiletto.
5. Insert an anterior chamber (AC) maintainer with balanced salt solution (BSS) or air.
6. Create inferior iridectomy.
7. Create a descemetorhexis (7-7.5 mm) under BSS or air in the AC.
8. Inject Vision Blue dye to the AC to make sure no remnants of the Descemet layer linger.
Note: In case of severe bleeding or uncontrolled vitreous in the AC, consider cancelling EndoArt® implantation.
9. Examine the label on the unopened package of EndoArt® implant for proper specification and expiration date.
10. After opening the cardboard storage container, remove the pouch with a primary container, containing EndoArt® implant.
11. Inspect the pouch carefully for any signs of damage. The primary container is sterile only if the pouch is not torn and properly sealed.
12. **Important:** Handle the primary container aseptically after opening the pouch. Remove the sterilized primary container from the pouch. Transport the primary container to a sterile environment.
13. Check to ensure there are no leaks from a primary container. If there is leak, do not use the implant.
14. Open the primary container and take the EndoArt® implant's holder out.

15. Open the EndoArt® implant's holder and carefully remove the implant using a sterile, atraumatic instrument.
16. Inspect EndoArt® implant for damage or defects prior to use.
17. Before insertion, ensure the implant is in its right orientation: EndoArt® has a character "(F)" marked on the peripheral surface that should appear correctly (Fig. 1); otherwise, the implant should be overturned.
18. Insert EndoArt® implant into the eye's anterior chamber with a spatula, preferably when the AC maintainer is running
19. After insertion, ensure again the correct placement and orientation of the implant; otherwise, the implant should be overturned.
20. Position the implant with the air/ C3F8 gas 10% (recommended) or SF6 gas 20% bubble, according to center of pupil.
21. Suture primary entrance.
22. Fill the AC 80% in volume with gas (IOP should be around 25mmHg. If needed, inject BSS/gas to reach the requested IOP, make sure air bubble is no more than 80%).
23. Suture the implant (single suture, nonabsorbable, at exactly 12 O'clock perpendicular to the center, according to Fig. 2). The perpendicular orientation of the suture is important to prevent decentration. The suture will be removed not sooner than 12 weeks from obvious attachment. Recommended removing before 6 months post good attachment.
24. Massage the cornea gently for enhancing corneal implant attachment, from the suture to 6 o'clock
25. In case iridectomy is not potent, apply pupil dilation drops to prevent pupillary block
26. After the last adjustment, wait for 20 minutes on the operating table with the speculum allowing the cornea to dry
27. Consider applying a contact lens to help in epithelial healing and decrease pain
28. Cover the eye with ventilated eye shield for 14 days.

29. Leave the patient in a supine position, facing upwards, for at least 4 hours after the operation.
30. In case of partial or full device detachment, re-bubbling procedure may be performed using 10% C3F8 gas or 20% SF6 gas, as per the physician's discretion. The re-bubbling procedure should be performed in an operation room, after which the patient should be placed in a supine position, facing upwards, for at least 4 hours after the operation.


Fig 1. EndoArt® implant orientation


Fig 2. EndoArt® suturing

10. POST OPERATION MEDICAL TREATMENT

- Use steroids for 12-24 weeks
- Use antibiotics for at least 1 month post operation, and as long as using bandage contact lens
- After re-bubbling use antibiotics and steroids for at least two weeks
- To minimize the risk of infection, use a bandage contact lens only for the shortest duration necessary.
- Post-operative follow-up should be focused on OCT corneal thickness measurements, in addition to standard of care of EK surgery follow-up.

11. GUIDANCE TO PATIENT

- Guide the patient to avoid touching/rubbing the eye during the first three months after implantation.
- Guide the patient to refer to the physician in the event of any sudden pain, decrease in vision, redness or secretion from the eye following the implantation.
- Guide the patient to remain with transparent eye shield for 14 days post procedure/re-bubbling
- Guide the patient not to fly or increase geographical altitude post air injection as long as a gas bubble exists in the operated eye. This may lead to increased IOP, pain and pathological damage to the eye.
- Guide the patient to lay on his back (face upwards) as much as possible in the first week post- op/re-bubble.
- Guide the patient to avoid any heavy lifting or sports activities of any kind for at least 2 weeks after the surgery.
- Guide the patient not to swim for 4 weeks after surgery
- Guide the patient to avoid activities that put him/her at risk for infection (like gardening, cleaning out stalls, attics, etc.) for the first month after surgery.

12. POTENTIAL COMPLICATIONS AND PROPOSED SOLUTIONS

Complication	Implant Position	Proposed Solution
Cornea		
Corneal edema	Implant is fully detached	Re-bubbling/new implantation
	Implant is partially detached	Re-bubbling procedure
Corneal thinning, risk of perforation	Implant is attached	Follow corneal thickness with OCT regularly until a stable thickness is established and then every 4 months. If corneal thickness does not stabilize and continues to decrease, a close follow-up should be performed. In case of corneal thinning below 300 microns, especially in combination with epithelial melting, consider removing the implant. Close post-operative follow-up should be performed until stabilization for at least 6 months.
Perforation	Implant is attached	Surgical correction
Haze	NA	Consider implant removal, DSAEK/DMEK procedure

Complication	Implant Position	Proposed Solution
Cornea		
Implant Migration	Not centered	Implant re-attachment/ re-bubbling/new implantation
Eye infection	NA	Treat as per standard of care
Epithelial bullae	Implant is attached	Consider medication for IOP reduction
	Implant is partially detached	Re-bubbling
Glaucoma		
Pupillary block	NA	Perform inferior iridotomy and IOP reduction medications
Cornea		
High intraocular pressure	NA	Treat as per standard of care
Iris adhesions and elevated IOP	NA	Provide standard glaucoma treatment or adhesion treatment surgery
Hypotony	NA	Make sure the entrance ports are closed or check for ciliary dialysis
Dry Eye		
Epiphora	NA	Treat as per standard of care (dry eye treatment)

13. STERILIZATION AND STORAGE

The device is sterilized by moist heat.

The device is non-pyrogenic.

To provide maximum protection, store the device in its original, unopened package at room temperature (up to 55°C).

Avoid excessive heat or cold.

The implant must be used before the expiration date printed on the package label.

14. DISCLAIMER OF WARRANTY

Note: EndoArt® referred to as “the product,” has been manufactured under carefully controlled conditions. EyeYon Medical has no control over product use and usage conditions. EyeYon Medical disclaims all warranties and representations, both expressed and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. EyeYon Medical shall not be liable to any person or entity to any, incidental or consequential damages, including, without limitation, for any medical expenses in connection to the product, caused by any use, defect, correspondence with description, non-infringement and/or any malfunction of the product, whether a claim for such damages is based upon representation, warranty, contract, tort, operation of law, statutory and/or otherwise. No person has any authority to bind EyeYon Medical to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

15. LABEL GRAPHICAL SYMBOLS DEFINITION

Use by

Keep dry

Lot number

Keep away
from sunlight

Do not reuse

Upper limit of
temperature

Do not resterilize

Manufacturer

Catalog number

Date of manufacture

Consult instructions for
use or consult electronic
instructions for use

Non-pyrogenic

Do not use if
package is
damaged and
consult instructions
for use

Authorized
representative in the
European Community/
European Union

Importer

Method of sterilization
using moist heat

Medical device

Double sterile
barrier system

Caution: This device
should be sold by and
used on the order of a
licensed physician

16. REPORTING OF SERIOUS INCIDENT

Serious incident relating to EndoArt® implant should be reported to EyeYon Medical Ltd and to the competent authority of the Member State in which the user and/or patient is located.

17. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The Summary of Safety and Clinical Performance (SSCP) is available upon request. To obtain a copy, please contact us at:

E-mail: **info@eye-yon.com**

Manufacturer:

EyeYon Medical Ltd.

3 Golda Meir St. Ness Ziona, 7403648, Israel

Email: **info@eye-yon.com**

European Union Authorized Representative:

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MedNet EC-REP IIb GmbH